



Simplicity in
site preparation




Tapered Implant Surgical Kits and Components Assembly

QNTSK Premium Kit




The Tapered Implant Surgical Tray

Simplicity In Site Preparation



NOTE: Graphic Illustrates Component Location For QNTSK40 Premium Kit

Tapered Implant Dense Bone Tap Kit



Note: Tapping in dense bone is required. Exceeding an insertion torque of more than 90 Ncm may deform or strip the driver tip or the implant's internal hex and may possibly delay the surgical procedure. Please refer to the Surgical Manual.

QNTSK Surgical Tray Contents

QNTSK40 Premium Kit

QNTSK40U Upgrade Kit

QNTSK20 Starter Kit

	Description	Item #			
1	Tapered Implant Depth Indicator 3.25 mm(D) X 8.5 mm(L) 3 Ea.	NTDI3285	.	.	.
1	Tapered Implant Depth Indicator 3.25 mm(D) X 10 mm(L)	NTDI3210	.	.	.
1	Tapered Implant Depth Indicator 3.25 mm(D) X 11.5 mm(L)	NTDI3211	.	.	.
1	Tapered Implant Depth Indicator 3.25 mm(D) X 13 mm(L)	NTDI3213	.	.	.
1	Tapered Implant Depth Indicator 3.25 mm(D) X 15 mm(L)	NTDI3215	.	.	.
2	Tapered Implant Depth Indicator 4 mm(D) X 8.5 mm(L) 3 Ea.	NTDI485	.	.	.
2	Tapered Implant Depth Indicator 4 mm(D) X 10 mm(L)	NTDI410	.	.	.
2	Tapered Implant Depth Indicator 4 mm(D) X 11.5 mm(L)	NTDI411	.	.	.
2	Tapered Implant Depth Indicator 4 mm(D) X 13 mm(L)	NTDI413	.	.	.
2	Tapered Implant Depth Indicator 4 mm(D) X 15 mm(L)	NTDI415	.	.	.
3	Tapered Implant Depth Indicator 5 mm(D) X 8.5 mm(L) 3 Ea.	NTDI585	.	.	.
3	Tapered Implant Depth Indicator 5 mm(D) X 10 mm(L)	NTDI510	.	.	.
3	Tapered Implant Depth Indicator 5 mm(D) X 11.5 mm(L)	NTDI511	.	.	.
3	Tapered Implant Depth Indicator 5 mm(D) X 13 mm(L)	NTDI513	.	.	.
3	Tapered Implant Depth Indicator 5 mm(D) X 15 mm(L)	NTDI515	.	.	.
4	Tapered Implant Depth Indicator 6 mm(D) X 8.5 mm(L) 3 Ea.	NTDI685	.	.	.
4	Tapered Implant Depth Indicator 6 mm(D) X 10 mm(L)	NTDI610	.	.	.
4	Tapered Implant Depth Indicator 6 mm(D) X 11.5 mm(L)	NTDI611	.	.	.
4	Tapered Implant Depth Indicator 6 mm(D) X 13 mm(L)	NTDI613	.	.	.
4	Tapered Implant Depth Indicator 6 mm(D) X 15 mm(L)	NTDI615	.	.	.
5	ACT® Pointed Starter Drill	ACTPSD	.	.	.
6	ACT Reusable Twist Drill 2 mm(D) X 10 mm(L)	ACT2010	.	.	.
6	ACT Reusable Twist Drill 2 mm(D) X 15 mm(L)	ACT2015	.	.	.
6	ACT Reusable Twist Drill 2 mm(D) X 20 mm(L)	ACT2020	.	.	.
7	Certain® Countersink Drill, 4 mm(D)	ICD100	.	.	.
8	Quad Shaping Drill 3.25 mm(D) X 8.5 mm(L)	QSD3285	.	.	.
8	Quad Shaping Drill 3.25 mm(D) X 10 mm(L)	QSD3210	.	.	.
8	Quad Shaping Drill 3.25 mm(D) X 11.5 mm(L)	QSD3211	.	.	.
8	Quad Shaping Drill 3.25 mm(D) X 13 mm(L)	QSD3213	.	.	.
8	Quad Shaping Drill 3.25 mm(D) X 15 mm(L)	QSD3215	.	.	.
9	Quad Shaping Drill 4 mm(D) X 8.5 mm(L)	QSD485	.	.	.
9	Quad Shaping Drill 4 mm(D) X 10 mm(L)	QSD410	.	.	.
9	Quad Shaping Drill 4 mm(D) X 11.5 mm(L)	QSD411	.	.	.
9	Quad Shaping Drill 4 mm(D) X 13 mm(L)	QSD413	.	.	.
9	Quad Shaping Drill 4 mm(D) X 15 mm(L)	QSD415	.	.	.
10	Quad Shaping Drill 5 mm(D) X 8.5 mm(L)	QSD585	.	.	.
10	Quad Shaping Drill 5 mm(D) X 10 mm(L)	QSD510	.	.	.
10	Quad Shaping Drill 5 mm(D) X 11.5 mm(L)	QSD511	.	.	.
10	Quad Shaping Drill 5 mm(D) X 13 mm(L)	QSD513	.	.	.
10	Quad Shaping Drill 5 mm(D) X 15 mm(L)	QSD515	.	.	.
11	Quad Shaping Drill 6 mm(D) X 8.5 mm(L)	QSD685	.	.	.
11	Quad Shaping Drill 6 mm(D) X 10 mm(L)	QSD610	.	.	.
11	Quad Shaping Drill 6 mm(D) X 11.5 mm(L)	QSD611	.	.	.

QNTSK40 Premium Kit

QNTSK40U Upgrade Kit

QNTSK20 Starter Kit

	Description	Item #			
11	Quad Shaping Drill 6 mm(D) X 13 mm(L)	QSD613	.	.	.
11	Quad Shaping Drill 6 mm(D) X 15 mm(L)	QSD615	.	.	.
12	Narrow Posterior Hexed Driver - 17 mm(L)	PHD00N	.	.	.
13	Posterior Large Hexed Driver - 17 mm(L)	PHD02N	.	.	.
14	Narrow Right Angle Large Hexed Driver Tip - 24 mm(L)	RASH3N	.	.	.
15	Certain Universal Ratchet Ext. Driver 4.1, 5, 6 mm(D) - Short	IRE100U	.	.	.
16	Certain Universal Ratchet Ext. Driver 4.1, 5, 6 mm(D) - Long	IRE200U	.	.	.
17	Direction Indicator 2 mm(D) X 10 mm(L) 3 Ea.	DI100	.	.	.
18	ACT Drill Extension	ACTDE	.	.	.
19	Handpiece Connector	MDR10	.	.	.
20	Implant Mount 3.4 mm(D) X 15 mm(L)	MMC15	.	.	.
21	Internal Connection Universal Placement Driver Tip (Long)	IIPDTUL	.	.	.
22	Internal Connection Universal Placement Driver Tip (Short)	IIPDTUS	.	.	.
23	Tapered Implant Dense Bone Tap 3.25 mm(D) X 8.5 mm(L)	NTAP3285	.	.	.
23	Tapered Implant Dense Bone Tap 3.25 mm(D) X 10 mm(L)	NTAP3210	.	.	.
23	Tapered Implant Dense Bone Tap 3.25 mm(D) X 11.5 mm(L)	NTAP3211	.	.	.
23	Tapered Implant Dense Bone Tap 3.25 mm(D) X 13 mm(L)	NTAP3213	.	.	.
23	Tapered Implant Dense Bone Tap 3.25 mm(D) X 15 mm(L)	NTAP3215	.	.	.
24	Tapered Implant Dense Bone Tap 4 mm(D) X 8.5 mm(L)	NTAP485	.	.	.
24	Tapered Implant Dense Bone Tap 4 mm(D) X 10 mm(L)	NTAP410	.	.	.
24	Tapered Implant Dense Bone Tap 4 mm(D) X 11.5 mm(L)	NTAP411	.	.	.
24	Tapered Implant Dense Bone Tap 4 mm(D) X 13 mm(L)	NTAP413	.	.	.
24	Tapered Implant Dense Bone Tap 4 mm(D) X 15 mm(L)	NTAP415	.	.	.
25	Tapered Implant Dense Bone Tap 5 mm(D) X 8.5 mm(L)	NTAP585	.	.	.
25	Tapered Implant Dense Bone Tap 5 mm(D) X 10 mm(L)	NTAP510	.	.	.
25	Tapered Implant Dense Bone Tap 5 mm(D) X 11.5 mm(L)	NTAP511	.	.	.
25	Tapered Implant Dense Bone Tap 5 mm(D) X 13 mm(L)	NTAP513	.	.	.
25	Tapered Implant Dense Bone Tap 5 mm(D) X 15 mm(L)	NTAP515	.	.	.
26	Tapered Implant Dense Bone Tap 6 mm(D) X 8.5 mm(L)	NTAP685	.	.	.
26	Tapered Implant Dense Bone Tap 6 mm(D) X 10 mm(L)	NTAP610	.	.	.
26	Tapered Implant Dense Bone Tap 6 mm(D) X 11.5 mm(L)	NTAP611	.	.	.
26	Tapered Implant Dense Bone Tap 6 mm(D) X 13 mm(L)	NTAP613	.	.	.
26	Tapered Implant Dense Bone Tap 6 mm(D) X 15 mm(L)	NTAP615	.	.	.
	Open End Wrench*	CW100	.	.	.
	Certain Driver Replacement O-Rings (10-Pack)*	IRORDR	.	.	.
	Tapered Implant Dense Bone Tap Kit Tray*	NTAPKT	.	.	.
	Ratchet Extension Short 6 mm(L)*	RE100	.	.	.
	Ratchet Extension Long 15 (L)*	RE200	.	.	.
	Ratchet Wrench*	WR150	.	.	.
	Base, Molded Insert and Smoked Lid*	QNTSKT	.	.	.
	Implant Depth Gauge 20 mm(L)*	DPO20	.	.	.

* Items not pictured in this brochure.

For more information, visit [ZimVie.com](https://www.zimvie.com)

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272
dentalCS@ZimVie.com



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZVINST0149 REV A 10/23 ©2023 ZimVie. All rights reserved.

