



**Mergence™  
Anterior  
Preparation  
Instrumentation**



Versatility, Flexibility and Simplicity




**zimmer** | spine  
*Confidence in your hands®*

# Mergence Anterior Preparation Instrumentation

## Access and Disc Removal

Set: 07.01368.400



## Endplate Preparation

Set: 07.01369.400



## Access and Disc Removal



Osteotome  
12.75mm 07.01220.001  
6.5mm 07.01220.002



Syfert Rongeur, 8 x 16mm  
360mm 07.01221.001



Cobb Elevator  
25mm 07.01212.001  
12mm 07.01212.002



Pituitary Rongeur, IVD  
6 x 400mm, Str 07.01222.001  
4 x 400mm, Str 07.01222.002



Kerrison Rongeur  
6mm, 400mm 07.01223.001  
4mm, 400mm 07.01223.002



Nerve Probe, Rt Angle, Blunt  
380mm 07.01225.001



Penfield Elevator, Blunt  
4mm, 380mm 07.01226.001



Graft Tamp 07.01224.001



Mallet 07.01237.001



## Endplate Preparation



Cup Curette, Straight, Serrated  
9.5mm 07.01213.001  
6.5mm 07.01227.001



Rake Curette, Angled  
12.75mm 07.01218.002



Cup Curette, Up-Angled  
9.5mm 07.01214.001  
6.5mm 07.01214.002



Rasp  
12.75mm 07.01219.001  
6.5mm 07.01219.002



Cup Curette, Down-Angled  
9.5mm 07.01215.001  
6.5mm 07.01215.002



Ring Curette  
9.5mm 07.01216.001  
6.5mm 07.01216.002



Curette, Oval  
12.75mm 07.01217.001



Rake Curette, Straight  
12.75mm 07.01218.001





## Mergence™ Anterior Preparation Instrumentation

The *Mergence* Anterior Preparation Instrumentation was designed by spine surgeons to provide a complete set of tools for access, discectomy and endplate preparation from an anterior approach. It complements the Zimmer Spine implant-specific instrument sets to provide the total solution for anterior spinal surgery.

- The comprehensive instrument set provides multiple options to accommodate varying anatomy and surgeon preference.
- The handles are contoured for a two-handed grip and have a comfortable and durable silicone surface for enhanced ergonomics.
- The instrument shafts have an anti-glare coating for improved visualization.
- The versatile *Mergence* Lumbar Distractor is available for anterior, anterior/lateral and lateral procedures.

### Other Products Available from Zimmer Spine

- *Puros*®-A Anterior Allograft
- TM-400 VBR Device
- *Trinica*® Anterior Lumbar Plate System
- *Mergence* Lumbar Distractor



Zimmer Spine, Inc.  
7375 Bush Lake Road  
Minneapolis, MN 55439-2027  
U.S.A.

Phone 952.832.5600  
or 800.655.2614  
Fax 952.832.5620  
or 800.430.9110  
[www.zimmerspine.com](http://www.zimmerspine.com)

L1414 Rev. B ©2008 Zimmer Spine, Inc. (8515-1007-00)

Contact your Zimmer Spine representative or visit us at [www.zimmerspine.com](http://www.zimmerspine.com)

